Sil-Select Plus™

EN	FR	DE	ES	IT	PT	GR	NL
----	----	----	----	----	----	----	----

Document ID: FP09_I13_Plus_R01_E.5

Update: 01/12/2021

MATERIAL INCLUDED/	MATERIEEL INCLUS/
IM LIEFERUNGSMANG ENTHALTENE MATERIALEN/	MATERIALE INCLUIDO/
MATERIAL INCLUDED/	MATERIALE INCLUS/
MATERIALE INCLUIDO/	ΥΛΙΚΟ ΠΡΟΪΗΜΑΜΒΑΝΕΤΑ/
MATERIAAL AANWZIG	

Catalogue number	Referenční číslo
Bestelningsnummer	Bestellnummer
Número de catálogo	Numéro de catalogue
Κατάλογος αριθμός	Αριθμός καταλόγου
Catalogusnummer	Catalogusnummer

SP008	8 x 2.5ml	Sil-Select Plus™ Lower layer
	8 x 2.5ml	Sil-Select Plus™ Upper layer
	8 x 8ml	Sil-Select Plus™ Sperm Washing/Insemination medium
SP016	2 x 20ml	Sil-Select Plus™ Lower layer
	2 x 20ml	Sil-Select Plus™ Upper layer
	5 x 20ml	Sil-Select Plus™ Sperm Washing/Insemination medium
SP008PG	8 x 2.5ml	Sil-Select Plus™ Lower layer
	8 x 2.5ml	Sil-Select Plus™ Upper layer
	8 x 8ml	Sil-Select Plus™ Sperm Washing/Insemination medium with phenol red and gentamicin
SP016PG	2 x 20ml	Sil-Select Plus™ Lower layer
	2 x 20ml	Sil-Select Plus™ Upper layer
	5 x 20ml	Sil-Select Plus™ Sperm Washing/Insemination medium with phenol red and gentamicin
SP025LD	1 x 1.25ml	Sil-Select Plus™ Lower layer
SP025LP	1 x 1.25ml	Sil-Select Plus™ Upper layer
SP042LD	1 x 1.42ml	Sil-Select Plus™ Lower layer
SP042LP	1 x 1.42ml	Sil-Select Plus™ Upper layer
SP020LD	1 x 20ml	Sil-Select Plus™ Lower layer
SP020LP	1 x 20ml	Sil-Select Plus™ Upper layer
SP050LD	1 x 500ml	Sil-Select Plus™ Lower layer
SP050LP	1 x 500ml	Sil-Select Plus™ Upper layer
SP050PRLD	1 x 500ml	Sil-Select Plus™ Lower layer with phenol red and gentamicin
SP050PRLP	1 x 500ml	Sil-Select Plus™ Upper layer with phenol red and gentamicin
SP090_ID	1 x 100ml	Sil-Select Plus™ 80%
	1 x 500ml	Sil-Select Plus™ 80%
	1 x 50ml	Sil-Select Plus™ 80% with gentamicin
	1 x 100ml	Sil-Select Plus™ 80% with gentamicin

CALCULATIONS OF G-FORCES

The g-force of your centrifuge can be calculated using this formula:

OF	g = 1.118 x r x rpm²
	rpm = Square root (g / (1.118 x r))
	r = radius of centrifuge in mm
	rpm = rotations per minute / 1000

Example 1	Example 2
r = 100 mm	r = 100 mm
rpm = 1800 rotations per minute	g = 350g
g = 1.118 x 100 x 324 = 362g	rpm = 350 / (1.118 x 100) = 177
	rpm = 1770 rotations per minute

WARNINGS BEFORE USE

- Do not use the product if:

- It becomes discoloured (if medium contains phenol red), cloudy, or shows any evidence of microbial contamination;
- seal of the container is opened or defect when the product is delivered;
- expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Product which includes gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.
- Keep in its original packaging until the day of use

METHOD

Importantly, each laboratory should consult its own optimized and validated laboratory procedures.
For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

FRESH SEMEN SAMPLES

Mix the density gradient bottles by 5 bottle inversions before use.

- Bring all components of the system and samples to room temperature or to 37 °C.
- Transfer 2.5ml of Sil-Select Plus Upper layer into a sterile disposable centrifuge tube.
- Using a syringe with needle, place 1.25ml of Sil-Select Plus Lower layer under the Upper layer. Take care that the two layers are distinctly separated. This is done by placing the tip of the needle on the bottom of the test tube and slowly dispensing the Sil-Select Plus Lower layer. This two layer gradient is stable for up to two hours.
- Gently place up to 2.5ml of liquefied semen onto the Upper layer using a transfer pipette or syringe. Do not use a higher volume than the volume of the individual gradient layers or more than 10⁶ cells.
- Centrifuge for 15 to 18 minutes at 350g to 400g. When this centrifugation is completed you may not be able to visibly see a pellet. If so, it is essential to continue the procedure with a second centrifugation of 3 to 5 minutes.
- Remove supernatant down to the pellet.
- Using a syringe, add 2.3ml of sperm washing medium (e.g. FertiCult Flushing medium or Si-Select Plus Sperm Washing/Insemination medium, FertiPro NV) and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300g. Higher sperm concentration will require the maximum 10 minutes centrifugation to ensure a complete and thorough sperm wash.

- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate medium.

If samples do not liquefy and therefore do not pass through the layers, increasing the centrifugal force up to, but no more than, 500g will help to separate the sperm.

FROZEN SEMEN SAMPLES

Mix the density gradient bottles by 5 bottle inversions before use.

- Bring all components of the system and samples to room temperature or to 37 °C.
- Transfer 2ml of Sil-Select Plus Upper layer into a sterile disposable centrifuge tube.
- Using a syringe with needle, place 1ml of Sil-Select Plus Lower layer under the Upper layer. Take care that the two layers are distinctly separated. This is done by placing the tip of the needle on the bottom of the test tube and slowly dispensing the Sil-Select Plus Lower layer. This two-layer gradient is stable for up to two hours.
- Separately thaw the thawed semen sample onto the Upper layer using a transfer pipette or syringe (0.5ml maximum).
- Centrifuge for 15-20 minutes at 350g.
- Remove supernatant down to no less than the 0.5ml mark above the pellet.
- Using a syringe, add 2.3ml of sperm washing medium (e.g. FertiCult Flushing medium or Sil-Select Plus Sperm Washing/Insemination medium, FertiPro NV) and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300g.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate medium.

STORAGE/DISPOSAL INSTRUCTIONS

- Store products without gentamicin between 2.25 °C. Once opened: store between 2-8 °C.
- Store products with gentamicin between 2-8 °C.
- The products can be used up to 7 days after opening when sterile conditions are maintained and the products are stored at 2-8 °C.
- Keep away from sunlight
- The products are stable after transport (max. 5 days) at elevated temperatures (≤ 37 °C)
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

PRECAUTIONS

- Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-evaluated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven viruses transmitted with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP of Sil-Select Plus and Sil-Select Plus Sperm Washing/Insemination medium describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by scanning the following QR codes:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

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Sil-Select Plus Upper layer/Lower layer/80% media:

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Sil-Select Plus Sperm Washing/Insemination medium:

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