QUALITY ASSURANCE CERTIFICATE



PureSperm® Wash

for washing the sperm pellet recovered from a PureSperm® density gradient, for use in the swim-up method and for extending semen or sperm pellet for use in IUI.

Containing: 100 mL Catalogue no: PSW-100

Batch:Production date:Expiry date:100PSWZD0808 April, 201907 April, 2020

We certify that this product was manufactured according to set directives for medical devices. The quality of this particular batch was tested and evaluated objectively by stringent, quality control procedures.

Quality Control procedures:	Accepted values:	Test results:
Physical analyses:		
pH at room temperature in air	7.3-8.5	Comply
Osmolality (mOsm)	290-300	Comply
Sterility and toxin analyses:		
Microbiological growth control	No growth	Comply
Endotoxins, quantitative LAL assay	< 1.00	Comply
(EU/mL)		
Biological analyses:		
Human sperm survival after 18 hrs	>70 %	Comply

Date of release: 12 August 2019

This document is electronically generated and does not bear a signature.

Paul V. Holmes, MSc, PhD, DrMedScMarina Danilova,C.E.O.Quality Assurance Coordinator

Nidacon International AB Flöjelbergsgatan 16B, SE-431 37 Mölndal, Sweden Tel:+46-31-703 06 30 Fax:+46-31-40 54 15

E-mail: contact@nidacon.com Web: www.nidacon.com